




Economically test devices and components for burst strength and leakage over time with a single, easy-to-use, reliable tester.



AVID Innovation's compact benchtop BLT200 series Burst and Leak tester combines high pressure burst testing with static pressure leak monitoring to deliver a complete pressure test solution in a single piece of test equipment.




By simplifying the design of the BLT product family, AVID Innovation has lowered the total cost of ownership while raising the system's reliability.

The intuitive touch screen interface is operator friendly and features a password protected setup mode for quick test parameter changes by qualified personnel.

As with all AVID Innovation equipment, customization is possible to solve your particular application needs. Custom options include, but are not limited to, higher pressures, data capture and storage, PLC and/or SCADA system interfacing. Please contact AVID Innovation to discuss custom options in further detail.

- Test pressures from 0-500 PSI
- Performs both Burst and Leak tests
- Modes include: Ramp Burst, Step Burst, Leak, Manual
- Ramp rate and Pressure control
- Easy to use Touch Panel operation
- Simple design provides low cost, high reliability
- Custom options include higher pressures, data capture, PLC and/or SCADA interface

Dimensions (in):



AVID Innovation is a leading supplier of high quality, reliable and cost-effective equipment and machine solutions to the medical devices industry.

Visit our website to learn more about both our standard and customized automation solutions.

AVID Innovation

Specifications:

Operating Pressure Range	0-500 PSI, Inert Gas Only
Pressure Accuracy	+/- .5% (2.5 PSI @ 500 PSI)
Pressure Resolution	0.1 psi
Pressure Fittings	1/4" NPT, inlet and outlet
Electrical	110VAC, 50/60 Hz, 1A
Dimensions	8" H x 12" W x 10" D
Weight	16 lb